ASMS 2019 ANNUAL CONFERENCE WORKSHOP



# MASS SPEC PACKAGES FOR R



## BEFORE YOU START, YOU NEED TO CONVERT YOUR DATA!

- Most mass spec instrument vendors save data files in proprietary formats that R can't directly read
- There are a variety of open formats for mass spec data,
   mzML is the main one currently in use
- Fortunately, open source tools exist that can convert most proprietary formats to mzML
- Unfortunately, this adds an extra processing step and mzML files can be large




### PROTEOWIZARD TO THE RESCUE

- Proteowizard is an open source project with lots of tools for mass spectrometry data processing and analysis <a href="http://proteowizard.sourceforge.net">http://proteowizard.sourceforge.net</a>
- Of particular note: MSConvert
   A GUI (and command line) tool for converting among
   MS data formats including proprietary → open formats
- Important: you have to use the Windows version to convert from proprietary to open formats




## THE MSCONVERT GUI IS EASY TO USE



There's also a command line version, but you must use Windows



# **EXAMPLE DATA CONVERSION WORKFLOW**






### BIOCONDUCTOR IS A GREAT RESOURCE FOR MS PACKAGES

- Bioconductor is a project and repository for 1700+ R packages for biological data analysis <a href="https://bioconductor.org">https://bioconductor.org</a>
- Includes lots of great packages for mass spectrometry data analysis
- Packages typically have reference documentation and examples (vignettes)




## **WORKING WITH BIOCONDUCTOR**



https://bioconductor.org



## **WORKING WITH BIOCONDUCTOR**



https://bioconductor.org



## INSTALLING BIOCONDUCTOR PACKAGES

- Installing packages from Bioconductor is different than the standard R method: don't use install.packages()
- First you need to install the Bioconductor framework;
  Within an R session, type:

```
if (!requireNamespace("BiocManager", quietly = TRUE))
  install.packages("BiocManager")
BiocManager::install()
```

you only need to do this once (or when you upgrade)

To install individual packages, type:

```
BiocManager::install("package_name")
```

Refer to "Install" on the Bioconductor site for more info



#### NOTABLE BIOCONDUCTOR MS PACKAGES

- MSnBase infrastructure for reading, processing and analyzing MS data
- mzR
  unified API for reading a variety of MS data formats
- MassSpecWavelet
  MS spectrum processing tools
- xcms comprehensive set of tools for MS analysis



### XCMS HAS A GREAT TUTORIAL ON LCMS PROCESSING/ANALYSIS

- LCMS data preprocessing and analysis with xcms <a href="https://bioconductor.org/packages/release/bioc/vignettes/xcms/inst/doc/xcms.html">https://bioconductor.org/packages/release/bioc/vignettes/xcms/inst/doc/xcms.html</a>
- Covers
  - Loading data ("on disk")
  - High-level data review
  - Chromatographic peak detection
  - Retention time alignment, and cross-experiment feature grouping



### **Demo Time**

Using example data from the OpenMS Tutorial data sets

https://www.openms.de/tutorials/

Download the Example Data

Example\_Data/Labelfree/datasets/lfq\_spikein\_dilution\_1.mzML